

NU SKIN DEVICES

BEAUTIFUL SKIN MAKES YOU FEEL UNSTOPPABLE.
WITH NU SKIN'S AT-HOME SKINCARE DEVICES, YOU HAVE THE FREEDOM
TO BUILD THE ROUTINE THAT HELPS YOU GET IT.



ageLOC LumiSpa iO & Activating Cleanser

Profound cleansing, exfoliating, and skin renewal

As a first step in your daily routine, morning and evening for two minutes each

The ultimate skin care optimisation. Visibly cleaner, softer, smoother, radiant and healthy-looking skin powered by smart personalisation, Intelligent Coaching and IoT connectivity to superpower your results.

Two-Sense Motion Technology, Smart IoT Technology



ageLOC Boost & Activating Serum

Brightening, glow

After cleansing and toning, once per day for two minutes.

Boosts your skin's natural beauty for brighter, plumper, and bouncier skin. The Activating Serum helps to improve the appearance of dark spots.

Variable Pulse Microcurrent Technology



ageLOC Galvanic Spa & Facial Gels

Anti-ageing skincare

After cleansing, 2-3 times per week

Firming and lifting your skin while removing impurities and helping to visibly reduce lines and wrinkles.

Direct/Galvanic Microcurrent



ageLOC WellSpa iO System

To revitalise, restore, and relax to help you look and feel your best.

Use three times per week to daily.

Support for visibly firmer and tighter skin for a more youthful appearance, overall feeling of wellbeing for better holistic wellness, and better pre and post workout feeling.

Personalised Adaptive Microcurrent Technology, Smart IoT Technology

WHAT IS IT FOR?

USAGE

BENEFITS

THE TECH BEHIND IT

Source Euromonitor International Limited; Retail Value RSP terms; all channels; 2017 to 2022. Beauty Systems are at-home Skin Care Beauty Devices that are exclusively paired or recommended to be used with a topical consumable of the same brand. Claim verification based on Euromonitor custom research and methodology conducted April - June of 2023. Sales of at-home skin care beauty devices includes sales of electric facial cleansers as defined in Passport database. This category does not include hair care/removal appliances, body shavers, and oral care appliances.



ageLOC LumiSpa iO
ACTIVATING CLEANSER

*INPUT SMART SKINCARE
OUTPUT RADIANT, HEALTHY-LOOKING SKIN*

WHAT MAKES AGELOC LUMISPA iO A SKINCARE SUPERPOWER?

ageLOC LumiSpa iO is a multi-tasking beauty solution. It features a proprietary Two-Sense Motion Technology that when used along with the ageLOC LumiSpa Activating Cleansers, thoroughly cleanses, gently exfoliates and delivers 7 amazing skin benefits. But that's not it – ageLOC LumiSpa iO connects to the Nu Skin Vera app to optimise your skin care journey with Customised Skin Care Routines, Intelligent Coaching, Awards, Selfie Timeline, and more!

BENEFITS

100% OF USERS HAD SMOOTHER, SOFTER AND VISIBLY CLEANER SKIN*1

53% TO 62% OF USERS HAD LESS VISIBLE PORES***1

62% MORE SKIN RADIANCE**1

71% REDUCTION IN APPEARANCE OF BLEMISHES²

CLARITY
FRESHNESS

ACTIVATING CLEANSERS

Specifically formulated (and tested) to be paired with your LumiSpa iO, the five ageLOC LumiSpa Activating Cleansers are made with different skin types in mind: Dry, Normal to Combination, Oily, Sensitive, and Blemish Prone. They enable precise cushioning, cleansing, and interaction with the skin, optimising the effectiveness of LumiSpa iO's proprietary Two-Sense Motion action.



HOW TO USE

To get started, pair your LumiSpa iO device with the Nu Skin Vera app on your phone. Once paired, select a pre-loaded routine or create your own on the app!

STEP 1

Wet your face with lukewarm water. Dispense an ample amount (bigger than an almond) of the LumiSpa Activating Cleanser that best suits your skin type onto your hand. Apply evenly to your face. Wet the Silicone Head and turn on your device.

STEP 2

Cleanse with your LumiSpa iO device for two minutes. Glide it slowly across your face in broad strokes. Don't scrub or press too hard. The device will vibrate if you are misusing it. It will pause every 30 seconds. This is your cue to move to another area of your face.

STEP 3

Rinse your face and gently pat it dry. Follow with next steps in your routine. For optimal results, use the ageLOC LumiSpa iO System for two minutes, twice per day.

STEP 4

Remove the Silicone Head from the device. Rinse it with water while rubbing to remove the residual cleanser. Dry thoroughly. Rinse device under water. Wipe the device dry.

STEP 5



SMART IoT FEATURES AND EYE CARE

Thanks to Bluetooth, the device can connect with the Nu Skin Vera app on your phone unlocking features like Intelligent Coaching, Customised Skin Care Routines, Selfie Timeline, Easy Reordering, usage insights, personalised product recommendations, awards, and more!

Complete your routine by caring for the delicate skin beneath and around your eyes. Use ageLOC LumiSpa iO Accent & IdealEyes with your device to exfoliate and invigorate your eye area.

¹Clinical studies conducted on ageLOC LumiSpa System. The results and efficacy can be extrapolated to ageLOC LumiSpa iO System.

^{**}/^{***}Results obtained from a third-party professional assessment based on a 12-week in-vivo consumer use study featuring 150 female subjects, from 35 to 65 years old, with different skin types. All subjects were instructed to use the LumiSpa System corresponding to their skin type, in place of their normal facial cleanser, morning and night, for 12 weeks. Results extracted above correspond to:

*The proportion of subjects who had visible skin improvement observed by the clinical grader at the end of the study (all combinations of skin types/LumiSpa Cleansers/Silicone Heads).

^{**}The observations made by the clinical grader after 12 weeks. The subjects were instructed to use the ageLOC LumiSpa with the Gentle Head and Sensitive Cleanser, Normal Head with Normal to Combo Cleanser and Firm Head with Oily Cleanser.

^{***}The observations made by the clinical grader, after 12 weeks on the facial skin after using one of the 3 following systems, Normal to Combo Skin/Normal Head (improvement of 53.6%), Dry Skin/Normal Head (improvement of 62.1%), Oily Skin/Firm Head (improvement of 56.5%), twice a day for 12 weeks.

²Study conducted on 30 subjects, males and females ages 18-45 with Blemish prone skin, using the Gentle Silicone Head and ageLOC LumiSpa Activating Cleanser - Blemish Prone for 12 weeks.



ageLOC Boost ACTIVATING SERUM

ENJOY BRIGHTER, PLUMPER SKIN IN JUST TWO MINUTES A DAY.

VARIABLE PULSE MICROCURRENT TECHNOLOGY

Intermittent, variable pulses are sent to the skin. Together with the Activating Serum, it optimises skin for a more youthful-looking complexion.



BENEFITS

HELPS EASE VISIBLE SIGNS OF PREMATURE SKIN AGEING.

BOOSTS YOUR SKIN'S NATURAL BEAUTY FOR BRIGHTER, PLUMPER, AND BOUNCIER LOOKING SKIN.

PROMOTES SKIN CELLULAR TURNOVER.

OFFERS GENTLE EXFOLIATION.

PATENT-PENDING VARIABLE PULSE MICROCURRENT TECHNOLOGY THAT OPTIMISES AND ACTIVATES SKIN FOR A MORE YOUTHFUL-LOOKING COMPLEXION.

THE ACTIVATING SERUM HELPS TO IMPROVE THE APPEARANCE OF DARK SPOTS.

ACTIVATING SERUM

Provides optimal conductivity for the Variable Pulse Microcurrent and skincare benefits from its key ingredients:

- ageLOC BLEND
- BRIGHTENING BLEND: SEA DAFFODIL EXTRACT(CEA), VITAMINE C DERIVATIVE, NIACINAMIDE
- CAVIAR LIME FRUIT EXTRACT
- PLUMPING BLEND: VARIUS FRUIT EXTRACTS, HYALURONIC ACID



HOW TO USE

After cleansing and toning, turn on device by swiping finger up full length of touch surface.

STEP 1

Dispense three pumps of ageLOC Boost Activating Serum onto the device's conductive surface. Spread over one half of the face.

STEP 2

Begin application by swiping finger up full length of touch surface. Glide conductive surface over face in an upward/outward motion.

STEP 3

Repeat Step 2 & Step 3 on other half of the face.

STEP 4

Massage excess Activating Serum into skin, then follow with your daily skincare regimen.

STEP 5

BEFORE/AFTER PICTURES¹



Before application

After 12 weeks



Before application

After 12 weeks

First visible results in just 7 days¹ and benefits that last throughout the day⁵.

¹ Results obtained from a third-party professional assessment based on a 12-week in vivo, single-blind consumer use study with 32 female subjects aged 25 to 40 who used the ageLOC Boost System once a day for 12 weeks.

⁵ Results obtained from a third-party professional assessment based on a 6-week in vivo, consumer use study with 20 female subjects, between 25 and 40 years old, who used ageLOC Boost Activating Serum together with the ageLOC Boost Device, 1/day, for 6 weeks.




ageLOC Galvanic Spa FACIAL GELS

LOOKING FOR THE SECRET TO YOUTHFUL-LOOKING SKIN? YOU'VE FOUND IT.

DIRECT/GALVANIC MICROCURRENT TECHNOLOGY

Like a set of magnets, the device uses direct currents to support the delivery of key ingredients to the skin. Using different settings, a negative current first moves the negatively charged Pre-Treat Gel to the skin, where it binds to impurities. A positive current then withdraws those impurities, while pushing the positively charged ageLOC ingredients in the Treatment Gel to the epidermis.



BENEFITS

HELP TO VISIBLY REDUCE LINES AND WRINKLES AND REVEAL A FIRMER, MORE LIFTED, CONTOURED, APPEARANCE.

TARGET THE SOURCES OF AGEING TO PRESERVE THE LOOK OF YOUTH AND REDUCE THE APPEARANCE OF AGEING.

REVIVE AND INVIGORATE TIRED, STRESSED SKIN.

HELPS COMBAT THE NEGATIVE EFFECTS OF OXIDATIVE STRESS CAUSED BY ENVIRONMENTAL STRESSORS.

GALVANIC CURRENT SUPPORTS THE SUPPLY OF BENEFICIAL INGREDIENTS TO THE EPIDERMIS.

AGELOC GALVANIC SPA FACIAL GELS

The secret to those amazing benefits? A combination of your device and these key ingredients:

- MARIGOLD AND SEA KELP EXTRACTS (IN THE PRE-TREAT GEL).
- ageLOC INGREDIENTS, ARGININE AND SKIN CONDITIONING AGENTS, GREEN TEA EXTRACT AND JANIA RUBENS EXTRACT (CEA) (IN THE TREATMENT GEL).



HOW TO USE

After cleansing, evenly apply the contents of one vial of Pre-Treat Gel onto the entire face and neck, avoiding eyes, eyelids, and lips.

STEP 1

Push the selection button until the number 1 is lit on the display. This application will take 2 minutes.

STEP 2

Thoroughly wet hands with NaPCA Moisture Mist or water. Hold the device with a damp hand while touching the chrome plate on the back and begin to slowly glide the Conductor head over your skin. Be sure to avoid sensitive areas like your eyes, eyelids, and lips.

STEP 3


When your galvanic application is over, the device will emit a long beep followed by two short beeps and automatically turn off. Remove any remaining gel from your skin and the device with a damp towel.

STEP 4

Repeat the entire process using a whole vial of ageLOC Treatment Gel and selecting application number 2 on your device. This application will take 3 minutes.

STEP 5

BEFORE/AFTER PICTURE⁶



Before application After **one** application

- 1 REDUCTION IN APPEARANCE OF CROW'S FEET.
- 2 REDUCTION IN APPEARANCE OF LINES AND WRINKLES IN THE UNDER-EYE AREA.
- 3 EYE BAG LIFT/IMPROVED SKIN FIRMNESS IN EYE AREA.
- 4 REDUCED APPEARANCE OF FINE LINES IN CHEEK AREA.
- 5 REDUCTION IN VISIBLE NASOLABIAL WRINKLES ('LAUGH LINES').
- 6 UPLIFT LOOK IN CHIN AREA/IMPROVED SKIN FIRMNESS.

After 4 weeks, study participants felt that the elasticity of their skin had improved by an average of 48%.⁶

⁶Results obtained from a 2020 third-party professional assessment based on a 4-week in vivo, single-blind, randomised consumer use study with 30 female subjects who used the system three times a week.



ageLOC WellSpa iO SYSTEM

TO REVITALISE, RESTORE, AND RELAX TO HELP YOU LOOK AND FEEL YOUR BEST.



HOW DOES AGELOC WELLSPA iO SYSTEM WORK?

ageLOC WellSpa iO is an exclusive beauty, wellness and self-care device system designed to help you look and feel your best. From the amazing ageLOC formulas to the unique design of the device, adaptive microcurrent and connected iO experience, the ageLOC WellSpa iO system provides revitalising, restoring and relaxing beauty and wellness benefits.



HOW TO USE

Revitalise routine

Apply a sufficient amount of ageLOC Body Serum onto the target areas of dry, clean.

Lightly glide the ageLOC WellSpa iO device over the skin up towards the centre of your body and back down in slow, circular movements.

Reapply as needed during the session and massage any remaining serum into the skin.

Restore routine

Apply a sufficient amount of ageLOC Body Activating Gel to the desired target areas on clean skin, and lightly glide on the body area in long, upward strokes.

Then rotate the ageLOC WellSpa iO device to move to the next section, until the entire body area has been covered.

Massage any remaining gel into the skin.

Relax routine

Apply a sufficient amount of ageLOC Body Activating Gel to the desired target areas on clean skin, and lightly glide the ageLOC WellSpa iO device on the targeted body areas in slow, upward circular strokes.

Rotate the ageLOC WellSpa iO device to move to the next section, until the entire body area has been covered.

Massage any remaining gel into the skin.

BENEFITS

REVITALISE ROUTINE

Designed for skin with an appearance of cellulite, it helps visibly firm and contour skin to deliver radiant, smoother looking skin.

Helps visibly firm and contour skin to deliver radiant, smoother looking skin.

RESTORE ROUTINE

Helps depuff the skin's appearance for a smooth look through its manual massage movement.

Supports overall feelings of wellbeing, helping you look better and feel more centred and balanced.

RELAX ROUTINE

Supports your body via a manual massage to help invigorate or relax sore muscles pre- or post-workout.

Helps you unwind and feel refreshed after a long day.

AGELOC BODY SERUM & AGELOC BODY ACTIVATING GEL

ageLOC Body Serum is applied during your Revitalise regimen to help visibly firm and tighten the skin, while its blend of hydrating ingredients support skin smoothness and a more youthful appearance.

ageLOC Body Activating Gel leaves a cooling sensation during your Restore and Relax ageLOC WellSpa iO routines, and helps target overall feelings of wellbeing and visible signs of skin ageing.



GOING FURTHER WITH THE FULL REVITALISE REGIMEN

For the full Revitalise Regimen, add ageLOC Body Polish and ageLOC Dermatic Effects. After cleansing with the ageLOC Body Polish in the shower and patting your skin dry, follow the usual steps of your Revitalise routine. Then finish with ageLOC Dermatic Effects using firm, upward circular motions.

SMART IOT FEATURES

Leads you every step of the way with real-time, guided video-tutorials.

Provides a personalised iO Score to help you track your usage and stay motivated to reach your beauty & wellness goals.

Offers Easy Reordering so you always have your products on hand.



BEFORE & AFTER



BEFORE & AFTER





Nu Skin Devices

COMPARISON TABLE

MAIN BENEFITS



AGELOC LUMISPA iO SYSTEM
(Device with ageLOC LumiSpa Activating Cleansers)



AGELOC LUMISPA iO ACCENT
(Device with ageLOC LumiSpa IdealEyes)



ageLOC BOOST SYSTEM
(Device with ageLOC Boost Activating Serum)



ageLOC GALVANIC FACIAL SPA SYSTEM
(Device with ageLOC Facial Gels)

UTILISES ADAPTIVE MICROCURRENT TECHNOLOGY				
WHEN USED WITH THE SPECIFICALLY FORMULATED TOPICAL, PROVIDES A COMFORTABLE MANUAL MASSAGE WHILE CONTRIBUTING TO SKIN STIMULATION				
GALVANIC CURRENTS SUPPORT THE SUPPLY OF KEY INGREDIENTS TO THE EPIDERMIS				X
WHEN USED WITH THE SPECIFICALLY FORMULATED TOPICAL, THE MICROCURRENTS OPTIMISE AND ACTIVATE SKIN			X	
TARGET VISIBLE SIGNS OF SKIN AGEING				X
TARGET THE PREMATURE SIGNS OF SKIN AGEING			X	
DECREASED APPEARANCE OF FINE LINES AND WRINKLES				X
MORE CONTOURED LOOKING SKIN				X
LIFTED UPPER EYELIDS		X		
FIRMER LOOKING SKIN		X (AROUND THE EYES)	X	X
IMPROVEMENT IN SKIN ELASTICITY				X
CONNECTED DEVICE WITH THE NU SKIN VERA APP	X	X		
CUSTOMISED ROUTINES	X	X		
GUIDANCE DURING USAGE	X	X		
PLUMPER, BOUNCIER LOOKING SKIN			X*	
RADIANT SKIN /BRIGHTNESS	X		X	
SKIN TONE EVENNESS			X	

*Linked to the sea daffodil extract, present in the Activating Serum.

**thanks to the movement of the device over the skin that helps massage the skin.

Disclaimer:

We have listed the main benefits, but please note that many of our devices + cosmetic products have secondary benefits as well.

MAIN BENEFITS



AGELOC LUMISPA iO SYSTEM
(Device with ageLOC LumiSpa Activating Cleansers)



AGELOC LUMISPA iO ACCENT
(Device with ageLOC LumiSpa IdealEyes)



ageLOC BOOST SYSTEM
(Device with ageLOC Boost Activating Serum)



ageLOC GALVANIC FACIAL SPA SYSTEM
(Device with ageLOC Facial Gels)

GENTLE EXFOLIATION TO PROMOTE SKIN CELLULAR TURNOVER	X	X	X	
SMOOTH & SOFT LOOKING SKIN	X	X	X	
THOROUGHLY CLEANSED SKIN	X			
IMPROVED PORE APPEARANCE	X		X	
REDUCED APPEARANCE OF EYE PUFFINESS/BAGS AND DARK CIRCLES		X		
APPEARANCE OF VOLUMINOUS, THICKER, HEALTHIER LOOKING HAIR				
STRENGTHENED AND FORTIFIED HAIR				
DESIGNED FOR SKIN WITH APPEARANCE OF CELLULITE				
IMPROVED APPEARANCE OF SAGGING SKIN				
TONED AND SCULPTED-LOOKING SKIN				
OVERALL FEELING OF WELLBEING				
IMPROVEMENT OF CIRCULATION IN THE SKIN LAYERS AND CONNECTIVE TISSUE**				
DEPUFFED SKIN'S APPEARANCE**				
SUPPORT THE BODY FOR A BETTER PRE AND POST WORKOUT FEELING				
HELP MUSCLES FEEL INVIGORATED AND READY TO GO				
HELP RELAX SORE MUSCLE				
OVERALL RELAXATION				

*Linked to the sea daffodil extract, present in the Activating Serum.

**thanks to the movement of the device over the skin that helps massage the skin.

Disclaimer:

We have listed the main benefits, but please note that many of our devices + cosmetic products have secondary benefits as well.

MAIN BENEFITS



ageLOC GALVANIC FACIAL SPA WITH TRU FACE LINE CORRECTOR



ageLOC GALVANIC FACIAL SPA WITH ageLOC NUTRIOL INTENSIVE SCALP & HAIR SERUM



ageLOC GALVANIC FACIAL SPA WITH ageLOC BODY SHAPING GEL



AGELOC WELLSPA iO - REVITALISE ROUTINE (DEVICE + AGELOC BODY SERUM)

UTILISES ADAPTIVE MICROCURRENT TECHNOLOGY				X
WHEN USED WITH THE SPECIFICALLY FORMULATED TOPICAL, PROVIDES A COMFORTABLE MANUAL MASSAGE WHILE CONTRIBUTING TO SKIN STIMULATION				X
GALVANIC CURRENTS SUPPORT THE SUPPLY OF KEY INGREDIENTS TO THE EPIDERMIS	X	X	X	
WHEN USED WITH THE SPECIFICALLY FORMULATED TOPICAL, THE MICROCURRENTS OPTIMISE AND ACTIVATE SKIN				
TARGET VISIBLE SIGNS OF SKIN AGEING	X	X (SCALP)	X	X
TARGET THE PREMATURE SIGNS OF SKIN AGEING				
DECREASED APPEARANCE OF FINE LINES AND WRINKLES	X (AROUND EYES/MOUTH/FOREHEAD)			
MORE CONTOURED LOOKING SKIN				X
LIFTED UPPER EYELIDS				
FIRMER LOOKING SKIN			X (THIGHS & ABDOMEN & UPPER ARMS)	X (THIGHS & ABDOMEN & ARMS)
IMPROVEMENT IN SKIN ELASTICITY				X
CONNECTED DEVICE WITH THE NU SKIN VERA APP				X
CUSTOMISED ROUTINES				X
GUIDANCE DURING USAGE				X
PLUMPER, BOUNCIER LOOKING SKIN				
RADIANT SKIN /BRIGHTNESS				X

*Linked to the sea daffodil extract, present in the Activating Serum.

**thanks to the movement of the device over the skin that helps massage the skin.

Disclaimer:

We have listed the main benefits, but please note that many of our devices + cosmetic products have secondary benefits as well.

MAIN BENEFITS



ageLOC GALVANIC FACIAL SPA WITH TRU FACE LINE CORRECTOR



ageLOC GALVANIC FACIAL SPA WITH ageLOC NUTRIOL INTENSIVE SCALP & HAIR SERUM



ageLOC GALVANIC FACIAL SPA WITH ageLOC BODY SHAPING GEL



AGELOC WELLSPA iO - REVITALISE ROUTINE (DEVICE + AGELOC BODY SERUM)

SKIN TONE EVENNESS				X
GENTLE EXFOLIATION TO PROMOTE SKIN CELLULAR TURNOVER				
SMOOTH & SOFT LOOKING SKIN	X			X
THOROUGHLY CLEANSED SKIN				
IMPROVED PORE APPEARANCE				
REDUCED APPEARANCE OF EYE PUFFINESS/BAGS AND DARK CIRCLES				
APPEARANCE OF VOLUMINOUS, THICKER, HEALTHIER LOOKING HAIR		X		
STRENGTHENED AND FORTIFIED HAIR		X		
DESIGNED FOR SKIN WITH APPEARANCE OF CELLULITE			X	X
IMPROVED APPEARANCE OF SAGGING SKIN				X
TONED AND SCULPTED-LOOKING SKIN				X
OVERALL FEELING OF WELLBEING				
IMPROVEMENT OF CIRCULATION IN THE SKIN LAYERS AND CONNECTIVE TISSUE**				
DEPUFFED SKIN'S APPEARANCE**				
SUPPORT THE BODY FOR A BETTER PRE AND POST WORKOUT FEELING				
HELP MUSCLES FEEL INVIGORATED AND READY TO GO				
HELP RELAX SORE MUSCLE				
OVERALL RELAXATION				

*Linked to the sea daffodil extract, present in the Activating Serum.

**thanks to the movement of the device over the skin that helps massage the skin.

Disclaimer:

We have listed the main benefits, but please note that many of our devices + cosmetic products have secondary benefits as well.

MAIN BENEFITS



AGELOC WELLSPA iO -
RESTORE ROUTINE
(DEVICE + AGELOC BODY
ACTIVATING GEL)



AGELOC WELLSPA iO -
RELAX ROUTINE
(DEVICE + AGELOC BODY
ACTIVATING GEL)

UTILISES ADAPTIVE MICROCURRENT TECHNOLOGY	X	X
WHEN USED WITH THE SPECIFICALLY FORMULATED TOPICAL, PROVIDES A COMFORTABLE MANUAL MASSAGE WHILE CONTRIBUTING TO SKIN STIMULATION	X	X
GALVANIC CURRENTS SUPPORT THE SUPPLY OF KEY INGREDIENTS TO THE EPIDERMIS		
WHEN USED WITH THE SPECIFICALLY FORMULATED TOPICAL, THE MICROCURRENTS OPTIMISE AND ACTIVATE SKIN		
TARGET VISIBLE SIGNS OF SKIN AGEING	X	X
TARGET THE PREMATURE SIGNS OF SKIN AGEING		
DECREASED APPEARANCE OF FINE LINES AND WRINKLES		
MORE CONTOURED LOOKING SKIN		
LIFTED UPPER EYELIDS		
FIRMER LOOKING SKIN		
IMPROVEMENT IN SKIN ELASTICITY		
CONNECTED DEVICE WITH THE NU SKIN VERA APP		
CUSTOMISED ROUTINES		
GUIDANCE DURING USAGE		
PLUMPER, BOUNCIER LOOKING SKIN		
RADIANT SKIN /BRIGHTNESS		
SKIN TONE EVENNESS		

*Linked to the sea daffodil extract, present in the Activating Serum.

**thanks to the movement of the device over the skin that helps massage the skin.

Disclaimer:

We have listed the main benefits, but please note that many of our devices + cosmetic products have secondary benefits as well.

MAIN BENEFITS



AGELOC WELLSPA iO -
RESTORE ROUTINE
(DEVICE + AGELOC BODY
ACTIVATING GEL)



AGELOC WELLSPA iO -
RELAX ROUTINE
(DEVICE + AGELOC BODY
ACTIVATING GEL)

GENTLE EXFOLIATION TO PROMOTE SKIN CELLULAR TURNOVER

SMOOTH & SOFT LOOKING SKIN

THOROUGHLY CLEANSED SKIN

IMPROVED PORE APPEARANCE

REDUCED APPEARANCE OF EYE PUFFINESS/BAGS AND DARK CIRCLES

APPEARANCE OF VOLUMINOUS, THICKER, HEALTHIER LOOKING HAIR

STRENGTHENED AND FORTIFIED HAIR

DESIGNED FOR SKIN WITH APPEARANCE OF CELLULITE

IMPROVED APPEARANCE OF SAGGING SKIN

TONED AND SCULPTED-LOOKING SKIN

OVERALL FEELING OF WELLBEING

IMPROVEMENT OF CIRCULATION IN THE SKIN LAYERS
AND CONNECTIVE TISSUE**

DEPUFFED SKIN'S APPEARANCE**

SUPPORT THE BODY FOR A BETTER PRE AND
POST WORKOUT FEELING

HELP MUSCLES FEEL INVIGORATED AND READY TO GO

HELP RELAX SORE MUSCLE

OVERALL RELAXATION

X

X

X

X

X

X

X

*Linked to the sea daffodil extract, present in the Activating Serum.

**thanks to the movement of the device over the skin that helps massage the skin.

Disclaimer:

We have listed the main benefits, but please note that many of our devices + cosmetic products have secondary benefits as well.